

---

## Auto-Mag® PCR-Pure

Cat # AMD-S002

Version 2.0

---

### Protocol for DNA double-sided size selection (96-well Plate Format)

This protocol can be used to generate DNA fragment libraries with a certain size range or to narrow fragment-size distribution. Through a two times selections process, the both smaller and larger fragments can be removed and the expected size range of DNA fragments can be obtained. This method is called double-sided size selection. The following protocol exemplifies a size selection of DNA fragment libraries with a size range of 250–400bp from a 50ul DNA samples. The information provided is an approximate guide. By altering the volume ratios of DNA fragment libraries and PCR-Pure reagents, other size ranges can be obtained. The optimal ratio for other size ranges should be determined empirically

#### Materials and Equipment to be supplied by User:

- 100 % Ethanol
- Elution Buffer (diH<sub>2</sub>O, 10mM Tris-HCl pH 8.0 or TE buffer)
- For 96 well format: 96 well cycling plate
- Appropriate magnetic separation device compatible with 96-well PCR plate
- Well calibrated pipettor and Disposable pipette tips

#### Before Starting

- Prepare fresh 80% ethanol (Prepare from absolute ethanol. Do not use denatured alcohol).
- Bring the Auto-Mag® PCR-Pure to room temperature for at least 30 min before use.
- Suspension Auto-Mag® PCR-Pure by vortexing.

#### Procedure

1. Completely float the Auto-Mag® PCR-Pure reagent until it appears homogeneous in color.
  2. For using the 96-well plate, a starting sample volume of 50µl is recommended. For smaller samples, add nuclease-free water to adjust the volume.
  3. For removing the large DNA fragments, first add 35µl of well dispensed Auto-Mag® PCR-Pure reagent into the wells containing 50µl DNA sample. Mix thoroughly by pipetting up and down 10 times or until homogenous.
- Note: The volume ratio of Auto-Mag® PCR-Pure reagent suspension to the sample is 0.7x.*
4. Incubate the mixture for 5 minutes at room temperature to allow the large DNA fragments are bound to Auto-Mag® PCR-Pure beads.
  5. Place the sample plate on the 96-well magnetic separation device and allow the magnetic beads to settle on the magnet for 5-10 minutes or until the supernatant turns completely clear.

---

6. Keep the sample plate on the magnet and transfer the ~85µl of supernatant into the well of a new plate and discard the beads that contain the unwanted large DNA fragments.

*Note: Do not disturb the attracted magnetic beads while aspirating the supernatant.*

7. Completely float the Auto-Mag® PCR-Pure reagent again until it appears homogeneous in color.

8. For the smaller DNA fragments cutoff, add 15µl of well dispersed Auto-Mag® PCR-Pure reagent into the well containing supernatants from step 6. Mix thoroughly by pipetting up and down 10 times or until homogenous.

*Note: The total volume ratio of Auto-Mag® PCR-Pure reagent suspension to the original sample is 1.0 x now; (35µl and 15µl to 50µl).*

9. Incubate the sample plate for 10 minutes at room temperature to allow the selected DNA fragment are bound to Auto-Mag® PCR-Pure beads.

10. Place the sample plate on the 96-well magnetic separation device and allow the magnetic beads to settle on the magnet for 5-10 minutes or until the supernatant turns completely clear. Carefully remove and discard the supernatant by pipetting.

11. Keep the sample plate on the magnet and add 200µl of 80 % ethanol to each well and incubate for 30 seconds at room temperature. Carefully remove and discard the supernatant by pipetting.

*Note: DNA Wash Buffer must be diluted with ethanol prior to use.*

12. Repeat Steps 11 for second 80% ethanol wash.

13. Keep the sample plate on the magnetic separation device and air dry the magnetic beads by incubating sample tube at room temperature for 5 minutes.

*Note: It is critical to completely remove all traces of liquid but take caution in not over drying the beads as this will reduce the yield.*

14. Remove the sample plate from the magnetic separation device. Add an appropriate volume (20-50µl) of Elution Buffer to each sample and mix by pipetting up and down 10 times. Incubate the sample plate for 5 minutes at room temperature.

15. Place the sample plate back on the magnetic separation device and wait 2 minutes or until the magnetic beads clear from solution.

16. Transfer all of supernatant to a new 96-well plate or an appropriate storage vessel. The ultra-pure DNA fragments (250-400bp) are ready for next step of NGS library preparation process or store at -20°C.