

**Product Information**

<b>Product Name:</b> Auto-Mag® Viral RNA/DNA Isolation Kit	
<b>Product No.</b> AMD-DR01	<b>Lot No.</b> DR0145216
<b>Manufacture Date:</b> 2023-10-17	<b>Expiry Date:</b> 2025-10-17
<b>Note:</b> <i>For Research Use Only. Not for use in diagnostic procedures</i>	

**Components in the package:**

Description	Storage Condition	Additional Infor.
Auto-Mag® C-5	2-8 °C	
VDR Buffer	15-25 °C	
DRW1 Buffer*	15-25 °C	
DRW2 Buffer*	15-25 °C	
RNA Elution Buffer	15-25 °C	
Proteinase K	2-8 °C	

**Quality Control:**

Auto-Mag® Viral RNA/DNA Isolation Kit is qualified by isolating viral RNA from 200ul serum of the spiked Human Coronavirus 229E following the protocols outlined in the manual.

Methods/ Specification:	Results
Purified RNA from Human Coronavirus 229E is detected by One-step RT PCR using viral gene specific primers and 40 cycles.	confirm
The CT value of positive sample must be < 35.	pass

**Quality authorized by QC:**

Signature 

**Date:** 2023-11-2

**ISO Certification:**

Manufactured by AMD Biotech Inc. in compliance with ISO 9001-2015 certified quality management system.

1860 Montreal Rd. Tucker, GA. 30084

US/Canada 1-404-290-5063 || Web. [www.amdbiotech.com](http://www.amdbiotech.com) || E-mail. support@amdbiotech.com