

Product Information

Product Name: Auto-Mag® DTR	
Product No. AMD-S001	Lot No. S00144568
Manufacture Date: 2022-01-17	Expiry Date: 2024-01-17
<i>For Research Use Only. Not for use in diagnostic procedures</i>	

Components in the package:

Subcomponents Description	Storage Condition	Additional Infor.
Auto-Mag® DTR	2-8 °C, (Do not freeze)	

Quality Control:

Auto-Mag® DTR has been tested by using A DNA sequencing reaction of a standard plasmid pUC19 following the protocols outlined in the manual and purified sequence products are analyzed on an ABI 3730 Sequencer.

Methods/ Specification:	Results
Conductivity (µS/cm) 106 µS/cm	pass
Average length of read	pass
No fluorescent contaminants	pass

Quality authorized by QC:

Signature 

Date: 2022-01-20

ISO Certification:

Manufactured by AMD Biotech Inc. in compliance with ISO 9001-2015 certified quality management system.