

**Product Information**

<b>Product Name:</b> Auto-Mag® X-Pure Select	
<b>Product No.</b> AMD-S003	<b>Lot No.</b> S00344852
<b>Manufacture Date:</b> 2022-10-18	<b>Expiry Date:</b> 2024-10-18
<i>For Research Use Only. Not for use in diagnostic procedures</i>	

**Components in the package:**

Subcomponents Description	Storage Condition	Additional Infor.
Auto-Mag® X-Pure Select	2-8 °C, (Do not freeze)	

**Quality Control:**

Auto-Mag® X-Pure Select has been tested by using the procedure of PCR reaction cleanup, and DNA fragment size-selection following the protocols outlined in the manual.

Methods/ Specification:	Results
The quality of the purified PCR is evaluated by agarose gel electrophoresis.	confirm
The yield and quality of DNA size selection is determined by TapeStation system	confirm

**Quality authorized by QC:**

Signature 

**Date:** 2022-10-25

**ISO Certification:**

Manufactured by AMD Biotech Inc. in compliance with ISO 9001-2015 certified quality management system.