

Product Information

Product Name: Auto-Mag® Tissue Total RNA Isolation Kit	
Product No. AMD-R023	Lot No. R02344747
Manufacture Date: 2022-07-05	Expiry Date: 2024-07-04
Note: <i>For Research Use Only. Not for use in diagnostic procedures</i>	

Components in the package:

Description	Storage Condition	Additional Infor.
Auto-Mag® C-5	2-8 °C	
RLB Buffer	15-25 °C	
DCE Buffer	15-25 °C	
MRW Buffer	15-25 °C	
RW2 Buffer	15-25 °C	
DNase I Digestion Buffer	15-25 °C	
DNase I	-20 °C	
Proteinase K	2-8 °C	
RNA Elution Buffer	15-25 °C	

Quality Control:

Auto-Mag® Tissue Total RNA Isolation Kit is qualified by isolating total RNA from 1×10^6 culture cells following the protocols outlined in the manual.

Methods/ Specification:	Results
A260/280 ratio of Purified total RNA must be between 1.8-2.2	pass
Total RNA is analyzed by gel electrophoresis. Two ribosomal RNA bands should be present, and the RNA should appear intact with no degradation.	confirm

Quality authorized by QC:

Signature 

Date: 2022-07-10

ISO Certification:

Manufactured by AMD Biotech Inc. in compliance with ISO 9001-2015 certified quality management system.