

Product Information

Product Name: Auto-Mag® Agarose Gel DNA Extraction Kit		
Product No. AMD-S010	Lot No. S01045021	
Manufacture Date: 2023-04-06	Expiry Date: 2025-04-06	
Note: For Research Use Only. Not for use in diagnostic procedures		

Components in the package:

Description	Storage Condition	Additional Infor.
Auto-Mag® D-1	2-8 °C	
GLB Buffer	15-25 °C	
GDW Buffer *	15-25 °C	
Elution Buffer	15-25 °C	

Quality Control:

Auto-Mag® Gel DNA Extraction Kit is qualified by extracting 0.5ug DNA ladder from agarose gel electrophoresis following the protocols outlined in the manual.

Methods/ Specification:	Results
The yield of the purified DNA was evaluated spectrophotometrically. DNA recovery must be ≥ 80%	pass
0.5 ug of extracted DNA ladder was compared with 0.5 µg reference DNA Ladder on agarose gel and all of fragment are purified.	confirm

Quality authorized by QC:

Signature

Date: 2023-04-06

ISO Certification:

Manufactured by AMD Biotech Inc. in compliance with ISO 9001-2015 certified quality management system.