

Product Information

Product Name: Auto-Mag® DNA Normalization Kit	
Product No. AMD-S006	Lot No. S00644761
Manufacture Date: 2022-07-19	Expiry Date: 2024-07-18
Note: <i>For Research Use Only. Not for use in diagnostic procedures</i>	

Components in the package:

Description	Storage Condition	Additional Infor.
Auto-Mag® C-7	2-8 °C	
NC Buffer	15-25 °C	
Auto-Mag® PCR-Pure	15-25 °C	
Elution Buffer	15-25 °C	

Quality Control:

Auto-Mag® DNA Normalization Kit is qualified by gDNA following the protocols outlined in the manual.

Methods/ Specification:	Results
Agarose gel electrophoresis of the normalized DNA must produce a band of the expected size, as compared to input DNA.	confirm
gDNA normalization should be recovery ~ 400ng of DNA with standard protocol.	confirm
Coefficient of Variation of normalization of recovered gDNA should be <10% between the 1500ng gDNA samples and the 6000ng samples	confirm

Quality authorized by QC:

Signature 

Date: 2022-07-31

ISO Certification:

Manufactured by AMD Biotech Inc. in compliance with ISO 9001-2015 certified quality management system.

1860 Montreal Rd. Tucker, GA. 30084

US/Canada 1-404-290-5063 || Web. www.amdbiotech.com || E-mail. support@amdbiotech.com