

Product Information

Product Name: Proteinase K, recombinant,	
Product No. AMD-B301	Lot No. B30145383
Manufacture Date: 04-01-2024	Expiry Date: 04-01-2027
Note: For Research Use Only. Not for use in diagnostic procedures	

Components in the package:

Description	Storage Condition	Additional Infor.
Proteinase K Solution (20mg/ml)	2-8 °C	

Quality Control:

Methods/ Specification/ Requirement	Results
Endodeoxyribonuclease. Incubation of supercoiled plasmid DNA with Proteinase K. Not detectable.	Confirm
Ribonucleases; Incubation of RNA transcript with Proteinase K. Not detectable.	confirm
Endonucleases; Incubation of DNA standard with Proteinase K. Not detectable.	confirm

Quality authorized by QC:
Signature
Date: 04-08-2024

ISO Certification:

Manufactured by AMD Biotech Inc. in compliance with ISO 9001-2015 certified quality management system.