

Product Information

Product Name: Linear Polyacrylamide (LPA)	
Product No. AMD-B315	Lot No. B31545338
Manufacture Date: 02-16-2024	Expiry Date: 02-06-2027
Note: For Research Use Only. Not for use in diagnostic procedures	

Components in the package:

Description	Storage Condition	Additional Infor.
Linear Polyacrylamide (LPA)	2-8 °C	

Quality Control:

Methods/ Specification/ Requirement	Results
Endonucleases activity: Sample incubation with plasmid DNA and analyzed by agarose gel electrophoresis. <10% change in intensity of full length, relative to controls.	Confirm
RNase activity: Sample incubation with RNA and analyzed by Tapestation, <10% change in intensity of full length, relative to controls.	confirm

Quality authorized by QC:
Signature

Date: 03-08-2024

ISO Certification:

Manufactured by AMD Biotech Inc. in compliance with ISO 9001-2015 certified quality management system.

1860 Montreal Rd. Tucker, GA. 30084

 US/Canada 1-404-290-5063 || Web. www.amdbiotech.com || E-mail. support@amdbiotech.com